

E-Z LabCheck

Features



Premium Kit Shown

- Base Kit comes in a rugged carrying case with one flow tube frame assembly, stand, flow meter support brackets, compression and hose barb fittings and conversion charts on CD;
- Up to three “E-Z Change” flow tubes can be used; allows quick interchangeability of flow tube assemblies to cover a wide range of applications;
- Rotating lens allows 180° view with magnification;
- Off-the-shelf reference and direct read metric scales available;
- Standard Kit with brass fittings and Buna-N seals, Premium Kit with stainless steel and Viton, Flow control valves available for each model;
- High accuracy calibrations:
 - Standard $\pm 3\%$ full scale
 - Optional $\pm .5\%$ of reading

APPLICATIONS

- Designed to enable users to comply with ISO 9000 self audits and company QA / QC requirements
- Secondary flow standard
- Ideal for bench testing and calibrating other flow measurement instruments
- Environmental/Analytical Instruments and Systems:
 - Verify sample acquisition rates
 - Set-up and balance flow circuits
 - Gas standard calibration accessory
 - Sample pump calibrations
- Instantaneously and accurately adjust flow rates
 - Diagnostics and troubleshooting



KEY INSTRUMENTS

250 Andrews Road • Trevoze, PA 19053 • 800-356-7483 • 215-357-0893 • Fax: 215-357-9239
www.keyinstruments.com • sales@keyinstruments.com

T-112

